Indirect Immunofluorescence Test for detection of
Taylorella equigenitalis

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Validated and certified by the OIE as fit for the purposes
defined in the kit insert. Registration number 20160111.

Pourquier* **IFI *Taylorella equigenitalis***

Test With Confidence™

 06-01410-07

IDEXX

Indirect Immunofluorescence Test for detection of *Taylorella equigenitalis*

For veterinary use only.

Name and Intended Use

Pourquier* IFI *Taylorella equigenitalis* is IDEXX's pool of monoclonal antibodies for the detection of *Taylorella equigenitalis* in genital swabs using immunofluorescence technique.

OIE Statement

The validation data for this kit have been certified by the OIE, based on expert review, as fit for the following purposes:

1. Certify freedom from infection or agent in individual animals or products for trade or movement purposes.
2. Estimate prevalence of infection to facilitate risk analysis (surveys, herd health schemes or disease control).
3. Control of infection in stallions and mares at the start of the breeding season.

Descriptions and Principles

Pourquier* IFI *Taylorella equigenitalis* uses the principle of immunofluorescence technique.

The swabs are fixed on slides with an acetone bath. The fixed swabs are placed in contact with the pool of monoclonal antibodies to allow bindings. A conjugate F(ab)'2 mouse anti-Immunoglobulin labelled with fluorescein isothiocyanate (FITC) is fixed on the monoclonal antibodies. The presence of *Taylorella equigenitalis* is demonstrated by observing the bacterial bodies presenting a typical fluorescence on cell wall, with a non fluorescent center. It shows no cross-reactivity with the following germs: *Actinobacillus equuli*, *Pseudomonas aeruginosa*, *Pasteurella multocida*, *Pasteurella haemolytica*, *Streptococcus equi*, *Staphylococcus aureus*, *Pseudomonas fluorescens* and *Klebsiella pneumonia*. It shows no cross-reactivity with *Tylorella asinigenitalis* which present antigenic community with *Taylorella equigenitalis*.

Reagents

Volume

1	Pool of monoclonal antibodies (3B6 1, 3B6 11, 7B7, 7C4)	1 x 1.2 mL
2	Conjugate (anti-mouse F(ab)'2 fragment affinity isolated FITC conjugate, diluted in Evans blue dye)	1 x 1.2 mL

Note: See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

Storage

Store all reagents away from light at $\leq -16^{\circ}\text{C}$ until the expiry date or up to 2 months at $2-8^{\circ}\text{C}$. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Phosphate buffer pH 7.2: the user can use any current formulation of phosphate buffer (without Ca^{2+} nor Mg^{2+}). The fluorescence of the antibody conjugated with the FITC is optimal when the pH is slightly alkaline, pH 7.2.
- Glycerol
- Acetone
- Precision Micropipettes and Multi-dispensing micropipettes
- Disposable pipette tips
- Humid chamber capable of maintaining a temperature of $+37^{\circ}\text{C}$ ($\pm 3^{\circ}\text{C}$)
- Dry incubator capable of maintaining a temperature of $+37^{\circ}\text{C}$ ($\pm 3^{\circ}\text{C}$)
- Alcohol 70%
- Fluorescent microscope with a 40x dry objective and an 100x oil-immersion objective (400-fold/1000-fold total magnification).
- Microscope slides with wells (2x5 wells) of 7mm diameter
- Cover-slides with dimensions adapted to microscope slides
- Staining glass basket with a panel for microscopic slides to achieve the various washes (ex. Dominique Dutscher with lid 105x85x75mm, ref. 068,506 and cart, ref. 068507)
- A McFarland scale
- Immersion oil for microscopy

Precautions and Warnings

- Handle all biological material as potentially infectious.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

Preparation of Reagents

Buffered glycerin, 9 volumes: Mix 1 volume of Phosphate buffer pH 7.2 with 9 volumes of glycerol.

Slide preparation for samples, positive and negative control

Use microscopic slides with wells cleaned with alcohol 70°C and dried.

In order to avoid cross-contamination between controls and field samples it is necessary to prepare the controls on separate slides.

For a positive control use any *Taylorella equigenitalis* reference strain or other strains isolated in the laboratory and correctly identified.

Use clean microscopic slides with wells, and add 30 μ l of a 0.5 McFarland suspension of *Taylorella equigenitalis* in the positive control well and 30 μ l of PBS in the negative control well.

Air dry the slides completely or incubate at 37°C (\pm 3°C) and fix them by immersion in acetone (under chemical hood or a glass box with cover) for 15 min at room temperature. Attention! Eliminate acetone after fixation.

Preparation of Samples

Genital swabs

Samples to be handled are in most cases swabs of the reproductive tract of stallions and mares (according to official regulations) placed in the transport medium with or without activated charcoal, such as Amies medium. The swabs should be suitably sealed with a stopper and send to the laboratory, in dark conditions within 3 days of sampling. Samples may be kept cool during transport.

Each swab is wiped by direct rotation on the slide wells (a complete rotation), being careful to avoid contamination between samples. Leave one empty well between sample wells for this purpose.

Dry the slides completely in a dry incubator at +37°C (\pm 3°C) for minimum 30 minutes.

Fix the slides by immersion in a pure acetone bath for 15 minutes at room temperature under a chemical hood. Attention! Eliminate acetone after the fixation of slides.

Note 1: If one sample contains a high number of *Taylorella equigenitalis*, some of these bacteria could be fixed on other wells of the same slide. For this reason, if other positive wells are observed on the same slide with a well containing a high number of *Taylorella*, this could be due to a cross-contamination between the samples.

Note 2: after fixation with acetone, the slides could be stored for few days in dark in a slide box avoiding the dust.

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Mix reagents by gentle inverting or swirling.

- 1 Dispense 30 μ l / spot of Pool of anti-Taylorella monoclonal antibodies.

- 2 Incubate for 30 minutes at +37°C (\pm 3°C) in a humid chamber.

- 3 Rinse the slides using a wash bottle with a trickle of phosphate buffer pH 7.2.

- 4 Immerse them in a phosphate buffer bath, agitate with a magnetic rod for 15 minutes.

- 5 Rinse the slides using a wash bottle with a trickle of distilled water.
Note: this last wash in distilled water is intended to prevent the formation of salt crystals resulting from the evaporation of the phosphate buffer.

- 6 Dry for 15 minutes in a dry incubator at +37°C (\pm 3°C).

- 7 Dispense 30 μ l / spot of Conjugate.

- 8 Incubate for 30 minutes at +37°C (\pm 3°C) in a humid chamber.

- 9 Rinse the slides using a wash bottle with a trickle of phosphate buffer pH 7.2.

- 10 Immerse them in a phosphate buffer bath, agitate with a magnetic rod for 15 minutes.

- 11 Rinse the slides using a wash bottle with a trickle of distilled water.
Note: this last wash in distilled water is intended to prevent the formation of salt crystals resulting from the evaporation of the phosphate buffer.

- 12 Dry for 15 minutes in a dry incubator at +37°C (\pm 3°C).
Note: the slide can be stored in the freezer at temperature \leq -16°C for later reading.

- 13 Set the slide with the buffered glycerin.

- 14** Read in a dark room using a fluorescent microscope with a total 400-fold and 1000-fold amplification. To focus the microscope bring the 40x objective as close as possible to the hydrophobic coating (without touching the area to be observed) using the rapid focusing wheel (to avoid crushing and breaking the blade). Place your eye on the eyepiece and slowly lower the deck with the rapid focusing wheel until the image becomes sharp (rough image). Improve the focus with the help of the fine focusing wheel. Read about 10 fields to verify fluorescent elements are present using the 40x dry objective. If fluorescent elements are present: confirm their nature with the 100x oil-immersion objective. If no fluorescent element is present: continue with the 40 dry objective and read the entire surface of the well by going back and forth from one edge of the well to the other and from top to bottom (or bottom to top) of the well.

Description of a *T. equigenitalis*-specific fluorescent element (typical fluorescence):

- Intense green fluorescence of the bacteria's cytoplasmic membrane with a non fluorescent center
- Coco-bacillary form which can also be more or less elongated.

Validation

The results of controls will validate the reaction as following:

- Negative control: no fluorescent bacteria should be observed
- Positive control: many bacterial bodies presenting a typical fluorescence

Interpretation of Results

The presence of *Taylorella equigenitalis* is demonstrated by observing at least one bacterial body presenting a typical fluorescence. The results are expressed by the "presence" or "absence" of specific fluorescent elements.

Summary of validation studies

Analytical characteristics

Repeatability: results comparable with the culture method for special agar with and without charcoal.

Analytical specificity: 100%

The tests are based on negative results of 18 pure cultures of bacterial strains other than *T. equigenitalis* which could be isolated from equine genital tract, 44 field isolates and one collection strain (CIP 107673T) of *T. asinigenitalis*.

Analytical sensitivity: 100%

The tests are based on the detection of all 706 pure cultures of field isolates of *T. equigenitalis* (506 field isolates from France and 200 field isolates from Holland) and 4 *T. equigenitalis* reference strains (NCTC 11225, ATCC 35865, ERC 7810381 and ERC 78107419).

The limit of detection (LOD) was evaluated in comparison with two culture methods (chocolate agar and selective Agar) in a mixture of *T. equigenitalis* and *Streptococcus zooepidemicus* in two different transport media (Amies with and without charcoal). The LOD of IFI was slightly superior in transport medium without charcoal (dilution 1/130,000 vs. 1/120,000) and in comparison with the culture method IFI was slightly superior to chocolate agar medium (dilution 1/130,000 vs. 1/100,000) and slightly inferior to selective agar culture medium (dilution 1/130,000 vs. 1/150,000).

Diagnostic characteristics

Diagnostic sensitivity (DSn) and specificity (DSp) estimates:

A. Field validation performed by OIE Reference Laboratory Netherlands (CVI Lelystad)

1) Validation based on individual horses

In the CVI validation study, from 730 horses tested, 12 horses were found positive by both culture and IDEXX IFI method. DSn and DSp of IDEXX IFI kit on individual horses considering culture as gold standard test:

<u>Test method under evaluation</u>	<u>IDEXX IFI</u>	<u>Target Species (Equine)</u>
Diagnostic sensitivity (Isolation and identification)	N DSn CI	12 (100%) 95% confidence interval: 73.5%-100%
Diagnostic specificity (Isolation and identification)	N DSp CI	718 (97.2%) 95% confidence interval: 95.7%-98.2%

2) Validation based on individual samples

A total of 2019 samples were tested by CVI and 1953 samples were found negative by both culture and IDEXX IFI test. The culture method yielded 19 positive samples and IFI method detected 65 samples as positive DSn and DSp of IDEXX IFI kit on individual samples considering culture as gold standard test:

<u>Test method under evaluation</u>	<u>IDEXX IFI</u>	<u>Target Species (Equine)</u>
Diagnostic sensitivity (Isolation and identification)	N DSn CI	19 (94.7%) 95% confidence interval: 73.9%-99.8%
Diagnostic specificity (Isolation and identification)	N DSp CI	2000 (97.6%) 95% confidence interval: 96.8%-98.2%

B. Field validation performed in case of non-experienced operator

In the field validation, from 2019 samples tested, 3 samples were found positive by culture and negative by IDEXX IFI test. These 3 samples were examined a second time by the same operators and in two slides specific *Taylorella equigenitalis* elements were found. The discordant results between the first reading and the re-examination of these 2 slides are explained by an error of non-experienced operators who did not detect the specific *Taylorella equigenitalis* elements at the routine examination of slides. Because of the fact that selective re-examination of samples (discrepant or discordance analysis) induces a bias in estimates of sensitivity and specificity, we calculate the DSn and DSp of the field validation after the first reading of slides. In this particular case, because of the low number of positive animals, the diagnostic sensitivity drops from 94% to 84% for individual samples and from 100% to 83% for individual horses.

1) Validation based on individual horses

DSn and DSp of IDEXX IFI kit on individual horses considering culture as gold standard test (non-experienced operator):

Test method under evaluation	IDEXX IFI	Target Species (Equine)
Diagnostic sensitivity (Isolation and identification)	N DSn CI	12 (100%) 95% confidence interval: 73.5%-100%
Diagnostic specificity (Isolation and identification)	N DSp CI	718 (97.2%) 95% confidence interval: 95.7%-98.2%

2) Validation on individual samples

DSn and DSp of IDEXX IFI kit on individual samples considering culture as gold standard test (non-experienced operator):

Test method under evaluation	IDEXX IFI	Target Species (Equine)
Diagnostic sensitivity (Isolation and identification)	N DSn CI	19 (84.2%) 95% confidence interval: 60.4%-96.6%
Diagnostic specificity (Isolation and identification)	N DSp CI	2000 (97.5%) 95% confidence interval: 96.7%-98.1%

Comparative performance:

Field validation performed by OIE Reference Laboratory Netherlands (CVI Lelystad) on individual samples.

All samples found positive by either culture or IFI were tested by PCR. From 22 individual samples found positive by PCR, 21 were positive by IFI DSn of IDEXX IFI on individual samples considering PCR as gold standard test.

Reproducibility

ANSES Dozulé (France) participates in some ring trials organized by one of the OIE laboratories for CEM (Animal Health Veterinary Laboratories Agency, GB). During the last session of June 2012 (report from the AHVLA dated July 6, 2012), the IDEXX IFI kit was used in parallel with the requested technique, bacteriology and PCR and ANSES obtained a 100% agreement with the requested technique.

The IDEXX kit was utilised for the ring trial in 2013 by 18 laboratories and in 2014 by 36 laboratories:

- In 2013 from 18 laboratories, 17 detected correctly all 6 positive swabs of the panel and one did not detect the highly diluted weak positive sample (100% correctly detected positive swabs).
- In 2014 from 36 laboratories, 34 detected correctly all 7 positive swabs of the panel and two laboratories did not detect the highly diluted weak positive sample (100% correctly detected positive swabs).

These results demonstrate a good reproducibility of the test to accurately detect the target organism even in a mixed bacterial flora which could be present in the reproductive tract of a mare and stallion.

For more details see the file "Ring trial for MCE by IFI: ANSES report 2013-2014."

Applications

2019 samples from 730 horses were tested with IDEXX IFI kit by CVI (Netherland) within the national screening program: export purposes (414 horses), screening (215 horses), certification artificial insemination service (91 horses) or for other reasons e.g. research/confirmation (10 horses).

In France, the IDEXX IFI kit has been used since 2013 principally for screening of stallions and mares at the start of breeding season (obligatory for “pure sang” and “trotteur français” races) and the movement of animals. All positive IFI animals should be confirmed by culture.

References

DA Gradinaru et al: Production and characterization of monoclonal antibodies against *Taylorella equigenitalis*, Vet Res (1997) 28, 65-76.

For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: idexx.com/contactlpd

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries

© 2018 IDEXX Laboratories, Inc. All rights reserved.

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IDEXX Montpellier SAS
IDEXX Europe B.V.