TeSeE™ Western Blot

Reagents for *in vitro* confirmation
of suspected TSE positive samples

Catalog # 355-1169

32 Tests

Validated and certified by the OIE for
the purposes defined in this insert

Registration number: 20090105

BIO-RAD

CONTENTS

- 1 - GENERAL INFORMATION
- 2 - ASSAY PRINCIPLE
- 3 - COMPOSITION OF THE KIT
- 4 - SAMPLES
- 5 - ASSAY PROCEDURE WITH MINI BLOT™ GEL
 - 5.1 Additional reagents and material required
 - 5.2 Preparation of reagents
 - 5.3 Sample purification
 - 5.4 Electrophoresis
 - 5.5 Protein transfer
 - 5.6 Immunoblotting
- 6 - ASSAY PROCEDURE WITH CRITERION™ XT GEL
 - 6.1 Additional reagents and material required
 - 6.2 Preparation of reagents
 - 6.3 Sample purification
 - 6.4 Electrophoresis
 - 6.5 Protein transfer
 - 6.6 Immunoblotting
- 7 - INTERPRETATION OF RESULTS
- 8 - PRECAUTIONS
- 9 - HYGIENE AND SAFETY INSTRUCTIONS
- 10 - REFERENCES

1 - GENERAL INFORMATION

Transmissible Spongiform Encephalopathies (TSE's) were first reported in the eighteenth century in sheep (Scrapie) and more recently in cervids such as deer and elk (Chronic Wasting disease, CWD) and cattle (Bovine Spongiform Encephalopathy, BSE). Humans are also susceptible to certain forms of TSE such as Kuru, Creutzfeldt-Jakob Disease (CJD) or Gerstmann-Sträussler-Scheinker Syndrome (GSS). The emergence of new variant Creutzfeldt-Jakob Disease (vCJD) in the human population has been strongly linked to the dietary intake of BSE-infected meat or meat products. One of the main characteristics of TSEs is a progressive accumulation in the central nervous system of an abnormal isoform of natural or cellular prion protein (PrP^c), termed PrP^{res}. This disease specific PrP^{res} is characterised by an increased resistance to proteases. The TeSeE™ WESTERN BLOT assay permits qualitative identification of PrP^{res} after proteolytic treatment which results in a reduced molecular weight fragment due to 'N' terminus truncation.

Active/passive surveillance programs have been conducted worldwide to detect BSE, scrapie or CWD in infected animals. Those programs have resulted in the identification of increased numbers of positive cases at the screening laboratories. Those positive samples (suspected animals) are then systematically confirmed as "TSE-infected" by the demonstration of typical spongiform changes with histopathology, or with the detection of abnormal PrP by Immunohistochemistry (IHC), or of Scrapie Associated Fibrils (SAFs) by electron microscopy. These above confirmation techniques require technical expertise for the interpretation of the results and are time consuming and expensive. Western blot technique can also be considered as an alternative method for confirmation of the TSE suspected samples.

The validation data for this kit have been certified by the OIE, based on expert review, as fit for the post-mortem detection of transmissible spongiform encephalopathies (TSEs) in cattle (bovine spongiform encephalopathy, BSE), in ovines and caprines (BSE and scrapie), and in cervids (Chronic Wasting Disease, CWD), and for the following purposes:

1. To confirm TSE suspected positive samples detected at the screening laboratories in countries with active/passive surveillance programmes. Any sample with a negative result according to the TeSeE™ WESTERN BLOT assay interpretation criteria, following a positive rapid test result, should be tested with one of the other OIE certified confirmatory methods, Immunohistochemistry (IHC) or SAF-Immunoblot;
2. To confirm the prevalence of infection with one of the TSE associated diseases (BSE, scrapie, CWD) in the context of an epidemiological survey in a low prevalence country;

3. To estimate prevalence of infection to facilitate risk analysis (e.g. surveys, implementation of disease control measures) and to assist the demonstration of the efficiency of eradication policies.

The TeSeE™ WESTERN BLOT assay is using the same assay principle as the Bio-Rad rapid assays (TeSeE™ SAP, TeSeE™ sheep/goat) that include the preliminary purification and concentration of the PrP^{res}, associated to a highly sensitive immunoblotting. Then, it can be used efficiently for confirmation of any TSE suspected samples and for typing of TSE strains in sheep.

2 - ASSAY PRINCIPLE

The TeSeE™ Western Blot assay allows the detection of PrP^{res} in nervous tissues (bovine, ovine, caprine, cervids, ...) or peripheral tissues (cervids) collected from infected animals.

The assay procedure begins with the digestion of cellular PrP protein (PrP^c), followed by purification and concentration of disease specific PrP^{res}. Detection of PrP^{res} is carried out by electrophoretic migration then immunoblotting using a monoclonal antibody highly specific for PrP^{res}.

The assay procedure includes the following steps:

- Sample homogenization,
- Digestion of PrP^c with proteinase K,
- Purification and concentration of PrP^{res},
- Electrophoresis and transfer onto a membrane,
- Immunoblotting.

3 - COMPOSITION OF THE KIT

Labelling	Type of reagents	Presentation	Storage
Grinding Tubes	Grinding tubes containing ceramic beads in a buffer solution ⁽¹⁾	1 bag (35 tubes)	+2°C to +25°C
A	Denaturing solution Ready to use	1 vial (20 ml)	+2°C to +25°C
B	Clarifying solution Colouring: bromophenol blue Ready to use	1 vial (20 ml)	+2°C to +25°C
PK	Proteinase K Colouring: phenol red	1 vial (0.5 ml)	+2°C to +8°C
Ab I	Primary antibody ⁽¹⁾ : anti-PrP monoclonal antibody (10x)	1 vial (8 ml)	+2°C to +8°C
Ab II	Secondary antibody ⁽¹⁾ : Sheep Anti-Mouse IgG-(H+L)-HRP (10x)	1 vial (10 ml)	+2°C to +8°C
Bl	Blocking solution ⁽¹⁾ (10x)	1 vial (10 ml)	+2°C to +8°C

⁽¹⁾ These reagents contain 0.1 % of ProClin® 300 (preservative).

4 - SAMPLES

The TeSeE™ WESTERN BLOT assay is suitable for the detection of TSEs in cattle (Bovine Spongiform Encephalopathy, BSE), in ovine and caprine (BSE and scrapie), and in cervids (Chronic Wasting Disease, CWD). This test can be processed directly from the same sample homogenate (grinding tube) prepared for Bio-Rad rapid testing (TeSeE SAP, TeSeE sheep/goat).

Bovine: purification of PrP^{res} is performed on samples from Central Nervous System (CNS). Since distribution of PrP^{res} is heterogeneous in central nervous system, obex area from brainstem must be preferably sampled for optimal detection.

Sampling syringe (Ref.: 355 1175) allows easy and rapid sampling of obex area in a secure way. Please refer to sampling protocol provided with sample syringes for detailed instructions on good sampling procedure.

Small ruminants: purification of PrP^{res} is performed on samples from Central Nervous System (CNS).

Samples are cut and weighed individually.

Cervids: purification of PrP^{res} is performed on samples from Central Nervous System (CNS) or peripheral tissues (lymphoid nodes). Samples are cut and weighed individually.

5 - ASSAY PROCEDURE WITH MINI BLOT™ GEL

5.1 - ADDITIONAL REAGENTS AND MATERIAL REQUIRED

5.1.1 - REAGENTS AND DISPOSABLES

Graduated pipettes (5, 10, 25 ml), conical tubes (50 ml), 2 ml polypropylene micro-test tubes with caps.

PARAFILM®M Sealing films.

Sample purification

Laemmli sample buffer	30 ml	Bio-Rad, cat. Nr. 161-0737
2-Mercaptoethanol	25 ml	Bio-Rad, cat. Nr. 161-0710
SDS	100 g	Bio-Rad, cat. Nr. 161-0301
Calibration syringes	200	Bio-Rad, cat. Nr. 355-1174

Gel electrophoresis

Acrylamide 40% 29:1	500 ml	Bio-Rad, cat. Nr. 161-0146
0.5 M Tris-HCl pH 6.8	1 L	Bio-Rad, cat. Nr. 161-0799
1.5 M Tris-HCl pH 8.8	1 L	Bio-Rad, cat. Nr. 161-0798
Bromophenol blue	10 g	Bio-Rad, cat. Nr. 161-0404
Sucrose	1 kg	Bio-Rad, cat. Nr. 161-0720
Ammonium persulfate	10 g	Bio-Rad, cat. Nr. 161-0700
TEMED	5 ml	Bio-Rad, cat. Nr. 161-0800
Tris/Glycine/SDS (running buffer) (10x)	1 L	Bio-Rad, cat. Nr. 161-0732
Kaleidoscope™ prestained standard	500 µl	Bio-Rad, cat. Nr. 161-0324
MagicMark™ XP Western Standard (Molecular weight standard)	250 µl	Invitrogen, cat. Nr. LC5602

Immunoblotting

Ethanol (Normapur)	1L	VWR, cat. Nr. 20821-296
Tris/CAPS (transfer buffer) (10x)	1L	Bio-Rad, cat. Nr. 161-0778
Filter paper (transfer paper for Mini Blot™ handcast gels)	50 sheets	Bio-Rad, cat. Nr. 170-3932
PVDF membrane (0.2 µm)	10 sheets	Bio-Rad, cat. Nr. 162-0175
Tween® 20	100 ml	Bio-Rad, cat. Nr. 170-6531
PBS (washing buffer) (10x)	1 L	Bio-Rad, cat. Nr. 161-0780
ECL (substrate for conjugate)	125 ml	Amersham, cat. Nr. RPN2109

ECL Hyperfilms (18 x 24 cm)	25 films	Amersham, cat. Nr. RPN2103K
Development folders	30 folders	Applied Biosystems, cat Nr. T2258
Kodak developing solution LX24	to 20 L	VWR or Kodak
Kodak fixative solution AL4	to 20 L	VWR or Kodak

5.1.2 - EQUIPMENT

Adjustable pipettes (10, 40, 200, 1000 µl),
 Graduated cylinder (1 L and 2 L), plastic forceps, trays, vortex.
 Exposure cassette and red light for film development.

Sample purification

TeSeE™ Precess 48™	Bio-Rad, cat. Nr. 359-0200
TeSeE™ Precess 24™	Bio-Rad, cat. Nr. 359-1070
Heater blocks (3 blocks)	Bio-Rad, cat. Nr. 358-9057
with Heating block – 20 tubes	Bio-Rad, cat. Nr. 358-9072
Centrifuge - 220/240 V	Bio-Rad, cat. Nr. 359-1396
Drum rotor	Bio-Rad, cat. Nr. 358-9189
Rotor adaptors - (x6)	Bio-Rad, cat. Nr. 358-9191

Gel electrophoresis

Mini-PROTEAN® Tetra Cell, electrophoresis module	Bio-Rad, cat. Nr. 165-8002
5 spacer plates	Bio-Rad, cat. Nr. 165-3312
PowerPac™ HC power supply: 100/120 V - 220/240 V	Bio-Rad, cat. Nr. 164-5052

Transfer

Trans-Blot® Cell	Bio-Rad, cat. Nr. 170-3946
------------------	----------------------------

Immunoblotting

Western Processor base	Bio-Rad, cat. Nr. 359-0098
Western Processor Mini Blot™ kit	Bio-Rad, cat. Nr. 170-3988

5.2 - PREPARATION OF REAGENTS

5.2.1 - SAMPLE PURIFICATION

- **Proteinase K**

Solution of proteinase K diluted in reagent A:

- ▶ 1 ml Reagent A
- ▶ 20 µl Proteinase K

Mix well by inverting until you obtain a homogeneous solution. After reconstitution, diluted proteinase K is stable 10 hours at room temperature (+18°C to +30°C).

- **Laemmli solution**

Solution of SDS + 2-Mercaptoethanol + Laemmli sample buffer:

- ▶ 0.6 g SDS
- ▶ 1.5 ml 2-Mercaptoethanol

Mix by inverting.

- ▶ 28.5 ml Laemmli sample buffer

Solution is aliquoted into 4 ml aliquots and stored at -20°C. Thawed aliquots can be re-frozen.

Note: It is recommended to prepare Laemmli solution one hour before use allowing SDS to be completely dissolved.

5.2.2 - ELECTROPHORESIS

- **Hand cast discontinuous acrylamide gel**

The gel must be 1.5 mm thickness.

Using the Mini Blot™ casting module, the resolving gel (13.5% acrylamide, pH 8.8) is cast first, once the resolving gel is polymerized the stacking gel is added (3% acrylamide, pH 6.8).

Resolving gel (1 gel)

- ▶ 2.8 ml Acrylamide 40%, 29:1
- ▶ 1.7 ml 1.5 M Tris-HCl buffer, pH 8.8 / SDS (1)
- ▶ 1.3 ml 50% sucrose solution (2)
- ▶ 2.5 ml distilled water

Mix by inverting.

- ▶ 43 µl 10% Ammonium persulfate (3)
- ▶ 9 µl TEMED

Pour 7 ml of the gel solution into the plates and retain the residual solution as a control of polymerization. Gently overlay to the top with 1 ml of 0.3 M Tris-HCl pH 8.8 / SDS buffer (4) so that the gel surface doesn't dry out. Let the gel polymerize for 15-20 minutes at room temperature (+18°C to +30°C). Check the residual solution is polymerized.

Invert the plate assembly to eliminate excess of buffer.

Stacking gel (1 gel)

- ▶ 4 ml 3% Acrylamide solution (7)
- ▶ 28 μ l 10% Ammonium persulfate (3)
- ▶ 6 μ l TEMED

Mix by inverting.

Gently pour the stacking gel onto the resolving gel and retain the residual solution as a control of polymerization. Position the comb, taking care not to trap any bubble in the well positions.

Let the gel polymerize for 5-10 minutes at room temperature (+18°C to +30°C). Check the residual solution is polymerized.

(1) Solution of 1.5 M Tris-HCl buffer, pH 8.8 / SDS

- ▶ 0.2 g SDS
- ▶ 50 ml 1.5 M Tris-HCl buffer pH 8.8

Solution can be stored at +2°C to +8°C for 2 weeks.

(2) Solution of 50% Sucrose

- ▶ 25 g Sucrose
- ▶ to 50 ml Distilled water

Sucrose solution can be stored at +2°C to +8°C for 2 weeks.

(3) Solution of 10% Ammonium persulfate

- ▶ 5 g Ammonium persulfate
- ▶ to 50 ml Distilled water

Ammonium persulfate solution is aliquoted and stored at -20°C. Thawed solution can be stored at +2°C to +8°C for 2 weeks.

(4) Solution of 0.3 M Tris-HCl buffer, pH 8.8 / SDS

- ▶ 40 ml Distilled water
- ▶ 10 ml 1.5 M Tris-HCl buffer pH 8.8 / SDS

Solution can be stored at +2°C to +8°C for 2 weeks.

(5) Solution of 0.5 M Tris-HCl buffer, pH 6.8 / SDS

- ▶ 0.2 g SDS
- ▶ 50 ml 0.5 M Tris-HCl buffer pH 6.8

Solution can be stored at +2°C to +8°C for 2 weeks.

(6) Solution of 1% Bromophenol Blue

- ▶ 0.5 g Bromophenol Blue
- ▶ 50 ml Distilled water

Bromophenol Blue solution can be stored at room temperature (+18°C to +30°C) for 6 months.

(7) Solution of 3% acrylamide

- ▶ 3.8 ml Acrylamide 40%, 29:1
- ▶ 10 ml 0.5 M Tris-HCl buffer pH 6.8 / SDS (5)
- ▶ 6 ml Sucrose 50% (2)
- ▶ 500 µl Bromophenol Blue 1% (6)
- ▶ to 50 ml Distilled water

Solution can be stored at +2°C to +8°C for 2 weeks.

• Kaleidoscope™ prestained standard

The Kaleidoscope™ prestained standard is prepared during the sample denaturation before loading on the acrylamide gel.

Prepare a 1/12 dilution in Laemmli [solution](#) (for example 10 µl of the Kaleidoscope™ prestained standard + 110 µl of Laemmli [solution](#)).

Please refer to the Kaleidoscope™ prestained standard insert for storage conditions.

• MagicMark™ XP Western Standard

The MagicMark™ XP molecular weight is prepared during the sample denaturation before loading on the acrylamide gel.

Prepare a 1/12 dilution in Laemmli [solution](#) (for example 10 µl of MagicMark™ XP + 110 µl of Laemmli [solution](#)).

Please refer to MagicMark™ XP insert for storage conditions.

• Mini Blot™ migration buffer

Solution of Tris-Glycine-SDS (1x).

Prepare a 1/10 dilution. **1 L of diluted buffer is required for 1 tank:**

- ▶ 900 ml Distilled water
- ▶ 100 ml Tris-Glycine-SDS buffer (10x)

Homogenize. Solution can not be stored.

5.2.3 - PROTEIN TRANSFER

- **Transfer buffer**

Solution of Tris/Caps-Ethanol 15%. **2.5 L is required for 1 transfer tank.**

- ▶ 750 ml Distilled water
- ▶ 150 ml Pure ethanol
- ▶ 100 ml Tris/CAPS (10x)

Homogenize. Solution can not be stored.

5.2.4 - Immunoblotting

- **Wash solution 1**

Solution of PBS (1x) + 0.1% Tween® 20. **Approximately 500 ml is required for the complete process of 1 membrane.**

- ▶ 900 ml Distilled water
- ▶ 100 ml PBS (10x)
- ▶ 1 ml Tween® 20

Thoroughly homogenize. Solution can be stored at +2°C to +8°C, overnight.

- **Wash solution 2**

Solution of PBS (1x). **Approximately 100 ml is required for the complete process of 1 membrane.**

- ▶ 900 ml Distilled water
- ▶ 100 ml PBS (10x)

Solution can be stored at room temperature (+18°C to +30°C) overnight.

- **Blocking solution**

During the transfer step, dilute the blocking solution (BI) 1/10 in Wash solution 1. **20 ml of diluted blocking solution (1x) is required for 1 membrane.**

- ▶ 18 ml Wash solution 1
- ▶ 2 ml Blocking solution (10x)

Homogenize by tube inverting.

- **Diluted primary antibody**

Just prior to use, dilute the primary antibody 1/10 in Wash solution 1.

15 ml of diluted primary antibody is required for 1 membrane.

- ▶ 13.5 ml Wash solution 1
- ▶ 1.5 ml Primary antibody (10x)

Homogenize by tube inverting.

- **Diluted secondary antibody (conjugate)**

Just before use, dilute the secondary antibody 1/10 in Wash solution 1.

20 ml of diluted conjugate is required for 1 membrane.

- ▶ 18 ml Wash solution 1
- ▶ 2 ml Secondary antibody (10x)

Homogenize by tube inverting.

- **ECL**

Substrate (ECL) must be prepared just before use. **1 ml of substrate is required for 1 membrane.**

- ▶ 0.5 ml Reagent 1
- ▶ 0.5 ml Reagent 2

Homogenize the solution.

- **Development solution**

- ▶ 800 ml Distilled water
- ▶ 200 ml Development product

Solution can be stored at room temperature (+18°C to +30°C), in a darkroom for 15 days maximum.

- **Fixative solution**

- ▶ 800 ml Distilled water
- ▶ 200 ml Fixative product

Solution can be stored at room temperature (+18°C to +30°C), in a darkroom for 15 days maximum.

5.3 - SAMPLE PURIFICATION

The TeSeE™ WESTERN BLOT assay can be processed directly from the same sample homogenate (grinding tube) prepared for Bio-Rad rapid tests (TeSeE SAP, TeSeE sheep/goat).

Sampling

For peripheral tissues (lymph nodes) insert one medium bead (Ref.: 355 1171) in the grinding tube, before to add the sample.

Take a mass of 350 mg ± 40 mg of nervous tissue (preferably in the obex area) or 200 mg ± 20 mg of peripheral tissue.

Deposit the sample in grinding tube, close firmly and proceed to the grinding step in the homogenizer (Ribolyser®, TeSeE™ PRECESS 48™ or TeSeE™ PRECESS 24™ - system).

Sample grinding

Place the tubes in the crown of the homogenizer.

Perform one agitation cycle with the following instrument parameters.

	Ribolyser®		TeSeE™ PRECESS 48™ or TeSeE™ PRECESS 24™	
	Nervous tissues	Peripheral tissues	Nervous tissues	Peripheral tissues
Time (sec.)	45	2 x 45 ⁽¹⁾	—	—
Speed	6.5	6.5	—	—
Program	—	—	Program 1	Program 2

When grinding is insufficient, another 1 or 2 agitation cycles can be performed⁽²⁾.

⁽¹⁾ ⁽²⁾ Wait a 5 minutes pause between the 2 agitation cycles.

Sample calibration

Remove the grinding tubes from the homogenizer, resuspend the homogenate by inverting before opening the tubes and aspirate 500 µl with the calibration syringe taking care to immerse the needle below the level of ceramic beads to avoid sampling tissue fragments.

Transfer each 500 µl sample into a 2 ml Eppendorf micro test-tube.

Note: At this stage, both grinding tubes after homogenisation and micro test-tubes after sample calibration can be stored, closed:

- At room temperature (+18°C to +30°C) for 15 hours.
- At +2°C to +8°C for 72 hours.
- At -20°C for 1 year. Frozen samples must be thawed at room temperature (+18°C to +30°C).

Samples can be submitted to a maximum of 3 freezing/thawing cycles.

Samples must always be homogenized by inverting before use.

Proteinase K treatment

Distribute 500 µl of reconstituted proteinase K solution (see paragraph 5.2.1) into each micro test-tube.

Homogenize the sealed tubes by inverting (10 times) and incubate at $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$ in a heating block incubator for 10 minutes.

Precipitation of PrPres with reagent B

Remove the tubes from the incubator. Open them and distribute 500 µl of reagent B into each tube. Homogenize by inverting until a homogeneous colour is obtained.

Concentration of the PrPres by centrifugation

Centrifuge the tubes for 7 minutes at 15000 g at 20°C .

Sample clarifying

Discard the supernatant by inverting over a waste container. Then dry the tubes by inverting onto absorbent paper for 5 minutes.

Distribute 100 µl of the Laemmli solution (see paragraph 5.2.1) into each micro test-tube.

Incubate for 5 minutes at room temperature ($+18^{\circ}\text{C}$ to $+30^{\circ}\text{C}$).

Completely resolubilize the pellet by aspiration/dispensing with a pipette.

Incubate for 5 minutes at $100^{\circ}\text{C} \pm 5^{\circ}\text{C}$ in a heating block incubator.

Remove the tubes from the incubator, homogenize by vortexing.

Centrifuge the tubes for 15 minutes at 15000 g at 20°C .

Transfer the supernatant to a new micro test-tube. Discard the tube containing the pellet.

At this stage, the supernatant can be stored frozen at -20°C for 24 hours; the samples must be thawed at room temperature ($+18^{\circ}\text{C}$ to $+30^{\circ}\text{C}$) prior to use.

5.4 - ELECTROPHORESIS

The TeSeE™ WESTERN BLOT assay can be used for both confirmation of TSE suspected samples or for strain typing in sheep.

The following procedure is applicable for confirmation of TSE suspected samples.

Please contact your Bio-Rad representative for instruction protocol in case of strain typing application.

Gel preparation

Place the acrylamide gels (see paragraph 5.2.2) in the migration tank. Pour the migration buffer (see paragraph 5.2.2) into the electrophoresis tank on each side of the gels, up to the top of the wells. Carefully remove the combs and rinse each well with migration buffer, using a pipette.

Sample loading

Heat the samples for 4 minutes at $100\text{ }^{\circ}\text{C} \pm 5^{\circ}\text{C}$ just prior to loading $15\text{ }\mu\text{l}$ /well. Load $15\text{ }\mu\text{l}$ of the diluted Kaleidoscope™ prestained standard and $15\text{ }\mu\text{l}$ of the diluted MagicMark™ XP (see paragraph 5.2.2).

Note: In case several gels are run at the same time, make sure that you stagger the loading of controls into different lanes for easy identification.

Differential migration of the samples

Run the gel at room temperature ($+18^{\circ}\text{C}$ to $+30^{\circ}\text{C}$) for 90 minutes at 150 V. The blue line must be out of the gel.

5.5 - PROTEIN TRANSFER

The transfer buffer must be prepared before the end of the sample migration (see paragraph 5.2.3).

Protein transfer preparation

Cut the membrane to the gel dimensions. Always use forceps when handling the membrane.

Immerse the membrane in pure ethanol for 15 seconds, rinse in distilled water for 5 minutes, then for 10 minutes in the transfer buffer.

Carefully remove the gel from the glass plates and let it equilibrate for 10 minutes in the transfer buffer.

Gel sandwich preparation

Soak filter paper and fibre pads in the transfer buffer.

Open the transfer cassette, with transparent side on the left. Respectively place on the transparent side a fiber pad, a filter paper, the membrane* and the gel*.

Complete with a filter paper then a fibre pad and close the cassette.

Immerse it in the transfer tank previously filled to the indicated limit with transfer buffer.

*Remove any air bubbles which may have formed.

Note: In case several membranes are processed at the same time, label each membrane in the corner.

Transfer onto the PVDF membrane

Agitate during the transfer by using a magnetic stirring bar and run for 60 minutes at 115 V.

5.6 - IMMUNOBLOTTING

a) Upon completion of the protein transfer, open the blotting assembly and remove the membrane for development. Quickly immerse the membrane in Wash solution 2 (see paragraph 5.2.4), then place it in ethanol for 10 seconds before rinsing for 5 minutes in distilled water.

Note: At this step, the membrane can be stored overnight in distilled water at +2°C to +8°C.

Let the membrane adjust to room temperature (+18°C to +30°C) before to start the immunoblotting.

b) Eliminate distilled water and incubate the membrane for 30 minutes in blocking solution (see paragraph 5.2.4). Incubate under medium agitation.

20 ml is sufficient for 1 membrane.

Note: from this step until the step g), the Bio-Rad Western Processor can be used for agitation and washing steps (refer to instruction manual for settings).

c) Eliminate the blocking solution and incubate the membrane in diluted **primary antibody** (see paragraph 5.2.4) for 30 minutes at room temperature (+18°C to +30°C) under medium agitation.

15 ml of diluted primary antibody is required for 1 membrane.

d) Eliminate the primary antibody solution and using Wash solution 1, briefly rinse the membrane, then wash twice for respectively 5 and 10 minutes, under fast agitation.

50 ml of Wash solution 1 is required for each cycle and for 1 membrane.

e) Eliminate Wash solution 1 and incubate the membrane for 20 minutes in diluted **secondary antibody** (see paragraph 5.2.4) at room temperature (+18°C to +30°C) under medium agitation.

20 ml of diluted secondary antibody is required for 1 membrane.

f) Eliminate the secondary antibody solution and using Wash solution 1, briefly rinse, then wash for respectively 5, 10 and 10 minutes under fast agitation.

50 ml of Wash solution 1 is required for each cycle and for 1 membrane.

g) Place the membrane in 50 ml of Wash solution 2 under slow agitation.

h) Drain the membrane on absorbent paper without blotting and place it in the plastic folder.

i) Add the ECL reagent (see paragraph 5.2.4). Eliminate the excess of reagent and air bubbles with absorbent paper. Place into the exposure cassette.

j) In a darkroom, cover the folder with a film and expose for 15 minutes. Film can be exposed longer or shorter time for optimal signal.

k) Immerse the film in developing solution for 45 seconds (see paragraph 5.2.4). Rinse in distilled water. Immerse the film in fixative solution until the film becomes completely transparent.

l) Wash with distilled water and let the film dry.

6 - ASSAY PROCEDURE WITH CRITERION™ XT GEL

6.1 - ADDITIONAL REAGENTS AND MATERIAL REQUIRED

6.1.1 - REAGENTS AND DISPOSABLES

Graduated pipettes (5, 10, 25 ml), conical tubes (50 ml), 2 ml polypropylene micro-test tubes with caps.

PARAFILM®M Sealing films.

Sample purification

Laemmli sample buffer	30 ml	Bio-Rad, cat. Nr. 161-0737
2-Mercaptoethanol	25 ml	Bio-Rad, cat. Nr. 161-0710
SDS	100 g	Bio-Rad, cat. Nr. 161-0301
Calibration syringes	200	Bio-Rad, cat. Nr. 355-1174

Gel electrophoresis

Criterion™ XT 12 % Bis-Tris	1 gel - 18 wells	Bio-Rad, cat. Nr. 345-0118
XT-MOPS (running buffer) (20 x)	500 ml	Bio-Rad, cat. Nr. 161-0788
Kaleidoscope™ prestained standard	500 µl	Bio-Rad, cat. Nr. 161-0324
MagicMark™ XP Western Standard (Molecular weight standard)	250 µl	Invitrogen, cat. Nr. LC5602

Immunoblotting

Ethanol (Normapur)	1L	VWR, cat. Nr. 20821-296
Tris/CAPS (transfer buffer) (10x)	1L	Bio-Rad, cat. Nr. 161-0778
Filter paper (transfer paper for Criterion™ XT precast gels)	50 sheets	Bio-Rad, cat. Nr. 170-4085
PVDF membrane (0.2 µm)	10 sheets	Bio-Rad, cat. Nr. 162-0175
Tween® 20	100 ml	Bio-Rad, cat. Nr. 170-6531
PBS (washing buffer) (10x)	1 L	Bio-Rad, cat. Nr. 161-0780
ECL (substrate for conjugate)	125 ml	Amersham, cat. Nr. RPN2109
ECL Hyperfilms (18 x 24 cm)	25 films	Amersham, cat. Nr. RPN2103K
Development folders	30 folders	Applied Biosystems, cat. Nr. T2258
Kodak developing solution LX24	to 20 L	VWR or Kodak
Kodak fixative solution AL4	to 20 L	VWR or Kodak

6.1.2 - EQUIPMENT

Adjustable pipettes (10, 40, 200, 1000 µl).

Graduated cylinder (1L and 2L).

Plastic forceps, trays, vortex.

Exposure cassette and red light for film development.

Sample purification

TeSeE™ Precess 48™,	Bio-Rad, cat. Nr. 359-0200
TeSeE™ Precess 24™,	Bio-Rad, cat. Nr. 359-1070
Heater blocks (3 blocks)	Bio-Rad, cat. Nr. 358-9057
with Heating block – 20 tubes	Bio-Rad, cat. Nr. 358-9072
Centrifuge - 220/240 V	Bio-Rad, cat. Nr. 359-1396
Drum rotor	Bio-Rad, cat. Nr. 358-9189
Rotor adaptors - (x6)	Bio-Rad, cat. Nr. 358-9191

Gel electrophoresis

Criterion™ XT Cell	Bio-Rad, cat. Nr. 165-6001
PowerPac™ HC power supply:	
100/120 V - 220/240 V	Bio-Rad, cat. Nr. 164-5052

Transfer

Criterion™ XT blotter	Bio-Rad, cat. Nr. 170-4070
-----------------------	----------------------------

Immunoblotting

Western Processor base	Bio-Rad, cat. Nr. 359-0098
Western Processor Criterion™ XT kit	Bio-Rad, cat. Nr. 170-3985

6.2 - PREPARATION OF REAGENTS

6.2.1 - SAMPLE PURIFICATION

- **Proteinase K**

Solution of proteinase K diluted in reagent A:

- ▶ 1 ml Reagent A
- ▶ 20 µl Proteinase K

Mix well by inverting until you obtain a homogeneous solution. After reconstitution, diluted proteinase K is stable 10 hours at room temperature (+18°C to +30°C).

- **Laemmli solution**

Solution of SDS + 2-Mercaptoethanol + Laemmli sample buffer:

- ▶ 0.6 g SDS
- ▶ 1.5 ml 2-Mercaptoethanol

Mix by inverting.

- ▶ 28.5 ml Laemmli sample buffer

Solution is aliquoted into 4 ml aliquots and stored at -20°C. Thawed aliquots can be re-frozen.

Note: It is recommended to prepare Laemmli solution one hour before use allowing SDS to be correctly dissolved.

6.2.2 - ELECTROPHORESIS

- **Kaleidoscope™ prestained standard**

The Kaleidoscope™ prestained standard is prepared during the sample denaturation before loading on the acrylamide gel.

Prepare a 1/12 dilution in Laemmli [solution](#), for example 10 µl of the Kaleidoscope™ prestained standard + 110 µl of Laemmli [solution](#).

Please refer to the Kaleidoscope™ prestained standard insert for storage conditions.

- **MagicMark™ XP Western Standard**

The MagicMark™ XP molecular weight is prepared during the sample denaturation before loading on the acrylamide gel.

Prepare a 1/12 dilution in Laemmli [solution](#), for example 10 µl of MagicMark™ XP + 110 µl of Laemmli [solution](#).

Please refer to MagicMark™ XP insert for storage conditions.

- **Criterion™ XT migration buffer**

Solution of MOPS (1x).

Prepare a 1/20 dilution. **1 L of diluted buffer is required for 1 tank:**

- ▶ 950 ml Distilled water
- ▶ 50 ml MOPS buffer (20x)

Homogenize. Solution can not be stored.

6.2.3 - PROTEIN TRANSFER

- **Transfer buffer**

Solution of Tris/CAPS-Ethanol 15%. **Approximately 2 L is required for 1 migration tank.**

- ▶ 750 ml Distilled water
- ▶ 150 ml Pure ethanol
- ▶ 100 ml Tris/CAPS (10x)

Homogenize. Solution can not be stored.

6.2.4 - IMMUNOBLOTTING

- **Wash solution 1**

Solution of PBS (1x) + 0.1% Tween® 20. **Approximately 1 L is required for the complete process of 1 membrane.**

- ▶ 900 ml Distilled water
- ▶ 100 ml PBS (10x)
- ▶ 1 ml Tween® 20

Thoroughly homogenize. Solution is stored at +2°C to +8°C, overnight.

- **Wash solution 2**

Solution of PBS (1x). **Approximately 200 ml is required for the complete process of 1 membrane.**

- ▶ 900 ml Distilled water
- ▶ 100 ml PBS (10x)

Solution is stored at room temperature (+18°C to +30°C) overnight.

- **Blocking solution**

During the transfer step, dilute the blocking solution (BI) 1/10 in Wash solution **1. 40 ml of diluted blocking solution is required for 1 membrane.**

- ▶ 36 ml Wash solution 1
- ▶ 4 ml Blocking solution (10x)

Homogenize by tube inverting.

- **Diluted primary antibody**

Just prior to use, dilute the primary antibody 1/10 in Wash solution 1.

30 ml of diluted antibody is required for 1 membrane.

- ▶ 27 ml Wash solution 1
- ▶ 3 ml Primary antibody (10x)

Homogenize by inverting.

- **Diluted secondary antibody (conjugate)**

Just prior to use, dilute the secondary antibody 1/10 in Wash solution 1.

40 ml of diluted conjugate is required for 1 membrane.

- ▶ 36 ml Wash solution 1
- ▶ 4 ml Secondary antibody (10x)

Homogenize by tube inverting.

- **ECL**

Substrate (ECL) must be prepared just prior to use. **2 ml of substrate is required for 1 membrane.**

- ▶ 1 ml Reagent 1
- ▶ 1 ml Reagent 2

Homogenize.

- **Development solution**

- ▶ 800 ml Distilled water
- ▶ 200 ml Development product

Solution is stored at room temperature (+18°C to +30°C), in a darkroom for 15 days maximum.

- **Fixative solution**

- ▶ 800 ml Distilled water
- ▶ 200 ml Fixative product

Solution is stored at room temperature (+18°C to +30°C), in a darkroom for 15 days maximum.

6.3 - SAMPLE PURIFICATION

The TeSeE™ WESTERN BLOT assay can be processed directly from the same sample homogenate (grinding tube) prepared for Bio-Rad rapid tests (TeSeE™ SAP, TeSeE™ sheep/goat).

Sampling

For peripheral tissues (lymph nodes) insert one medium bead (Ref.: 355 1171) in the grinding tube, before to add the sample.

Take a mass of 350 mg ± 40 mg of nervous tissue (preferably in the obex area) or 200 mg ± 20 mg of peripheral tissue.

Deposit the sample in grinding tube, close firmly and proceed to the grinding step in the homogenizer (Ribolyser®, TeSeE™ PRECESS 48™ or TeSeE™ PRECESS 24™ - system).

Sample grinding

Place the tubes in the crown of the homogenizer.

Perform one agitation cycle with the following instrument parameters.

	Ribolyser®		TeSeE™ PRECESS 48™ or TeSeE™ PRECESS 24™	
	Nervous tissues	Peripheral tissues	Nervous tissues	Peripheral tissues
Time (sec.)	45	2 x 45 ⁽¹⁾	—	—
Speed	6.5	6.5	—	—
Program	—	—	Program 1	Program 2

When grinding is insufficient, another 1 or 2 agitation cycles can be performed⁽²⁾.

⁽¹⁾ ⁽²⁾ Wait a 5 minutes pause between the 2 agitation cycles.

Sample calibration

Remove the grinding tubes from the homogenizer, resuspend the homogenate by inverting before opening the tubes and aspirate 500 µl with the calibration syringe taking care to immerse the needle below the level of ceramic beads to avoid sampling tissue fragments.

Transfer each 500 µl sample into 2 ml Eppendorf micro test-tube.

Note: at this stage, both grinding tubes after homogenisation and micro test-tubes after sample calibration can be stored, closed:

- At room temperature (+18°C to +30°C) for 15 hours.
- At +2°C to +8°C for 72 hours.
- At -20°C for 1 year. Frozen samples must be thawed at room temperature (+18°C to +30°C).

Samples can be submitted to a maximum of 3 freezing/thawing cycles.

Samples must always be homogenized by inverting before use.

Proteinase K Treatment

Distribute 500 µl of reconstituted proteinase K solution (see paragraph 6.2.1) into each micro test-tube.

Homogenize the sealed tubes by inverting (10 times) and incubate at 37°C ± 2°C in a heating block incubator for 10 minutes.

Precipitation of PrP^{res} with reagent B

Remove the tubes from the incubator. Open and distribute 500 µl of reagent B into each tube. Homogenize by inverting until a homogeneous colour is obtained.

Concentration of the PrP^{res} by centrifugation

Centrifuge the tubes for 7 minutes at 15000 g at 20°C.

Sample clarifying

Discard the supernatant by inverting over a waste container. Then dry the tubes by inverting onto absorbent paper for 5 minutes.

Distribute 100 µl of the Laemmli solution (see paragraph 6.2.1) into each micro test-tube.

Incubate for 5 minutes at room temperature (+18°C to +30°C).

Completely resolubilise the pellet by aspiration/dispensing with a pipette.

Incubate for 5 minutes at 100°C ± 5°C in a heating block incubator.

Remove the tubes from the incubator, homogenize by vortexing.

Centrifuge the tubes for 15 minutes at 15000 g at 20°C.

Transfer the supernatant to a new micro test-tube. Discard the tube containing the pellet.

At this stage, the supernatant can be stored frozen at -20°C for 24 hours; the samples must be thawed at room temperature (+18°C to +30°C) prior to use.

6.4 - ELECTROPHORESIS

The TeSeE™ WESTERN BLOT assay can be used for both confirmation of TSE suspected samples or for strain typing in sheep.

The following procedure is applicable for confirmation of TSE suspected samples.

Please contact your Bio-Rad representative for instruction protocol in case of strain typing application.

Gel preparation

Remove the plastic band on the bottom of the plastic plate and place the acrylamide gels (see paragraph 6.2.2) in the migration tank. Pour the migration buffer (see paragraph 6.2.2) on each side of the gel up to the top of the wells and into the electrophoresis tank. Carefully remove the combs and rinse each well with migration buffer, using a pipette.

Sample loading

Heat the samples for 4 minutes at $100\text{ }^{\circ}\text{C} \pm 5^{\circ}\text{C}$ just prior to loading $15\text{ }\mu\text{l}$ /well.

Load $15\text{ }\mu\text{l}$ of the diluted Kaleidoscope™ prestained standard and $15\text{ }\mu\text{l}$ of the diluted MagicMark™ XP (see paragraph 6.2.2).

Note: In case several gels are run at the same time, make sure that you stagger the loading of controls into different lanes for easy identification.

Differential migration of the samples

Run the gel at room temperature ($+18^{\circ}\text{C}$ to $+30^{\circ}\text{C}$) for 50 minutes at 200 V.

6.5 - PROTEIN TRANSFER

The transfer buffer must be prepared before the end of the sample migration (see paragraph 6.2.3).

Protein transfer preparation

Cut the membrane to the gel dimensions. Always use forceps when handling the membrane.

Immerse the membrane in pure ethanol for 15 seconds, rinse in distilled water for 5 minutes, then for 10 minutes in the transfer buffer.

Carefully remove the gel from the plastic plates and let it equilibrate for 10 minutes in the transfer buffer.

Gel sandwich preparation

Soak filter paper and fibre pads in the transfer buffer.

Open the transfer cassette, with red side on the left. Respectively place on the red side a fiber pad, a filter paper, the membrane* and the gel*.

Complete with a filter paper then a fibre pad and close the cassette.

Immerse it in the transfer tank, previously filled to the indicated limit with transfer buffer. A frozen ice pack is added prior to fill the tank.

*Remove any air bubbles which may have formed.

Note: In case several membranes are processed at the same time, label each membrane in the corner.

Transfer onto the PVDF membrane

Agitate during the transfer by using a magnetic stirring bar and run for 60 minutes at 115 V.

6.6 - IMMUNOBLOTTING

a) Upon completion of the protein transfer, open the blotting assembly and remove the membrane for development. Quickly immerse the membrane in Wash solution 2 (see paragraph 6.2.4), then place it in ethanol for 10 seconds before rinsing for 5 minutes in distilled water.

Note: At this step, the membrane can be stored overnight in distilled water at +2°C to +8°C.

Let the membrane adjust to room temperature (+18°C to +30°C) before to start the immunoblotting.

b) Eliminate distilled water and incubate the membrane for 30 minutes in blocking solution (see paragraph 6.2.4). Incubate under medium agitation.

40 ml is required for 1 membrane.

Note: from this step until step g), the Bio-Rad Western Processor can be used for agitation and washing steps (refer to instruction manual for settings).

c) Eliminate the blocking solution and incubate the membrane in diluted **primary antibody** (see paragraph 6.2.4) for 30 minutes at room temperature (+18°C to +30°C) under medium agitation.

30 ml of diluted primary antibody is required for 1 membrane.

d) Eliminate the primary antibody solution and using Wash solution 1, briefly rinse the membrane, then wash twice for respectively 5 and 10 minutes, under fast agitation.

100 ml of Wash solution 1 is required for each cycle and for 1 membrane.

e) Eliminate the Wash solution 1 and incubate the membrane for 20 minutes in diluted **secondary antibody** (see paragraph 6.2.4) at room temperature (+18°C to +30°C) under medium agitation.

40 ml of diluted secondary antibody is required for 1 membrane.

f) Eliminate the secondary antibody solution and using Wash solution 1, briefly rinse for respectively 5, 10 and 10 minutes under fast agitation.

100 ml of Wash solution 1 is required for each cycle and for 1 membrane.

g) Place the membrane in 100 ml of Wash solution 2 under slow agitation.

h) Drain the membrane on absorbent paper without blotting and place it in the plastic folder.

i) Add the ECL reagent (see paragraph 6.2.4). Eliminate the excess of reagent and air bubbles with absorbent paper. Place into the exposure cassette.

j) In a darkroom, cover the folder with a film and expose for 15 minutes. Film can be exposed longer or shorter time for optimal signal.


k) Immerse the film in developing solution for 45 seconds (see paragraph 6.2.4). Rinse in distilled water. Immerse the film in fixative solution until the film becomes completely transparent.

l) Wash with distilled water and let the film dry.

7 - INTERPRETATION OF RESULTS

Figure 1 shows the expected band patterns for TSE negative samples, TSE positive samples in various animal species and molecular weight controls (positions 3 and 9).

Figure 1

Negative samples (positions 1 and 10) were treated with proteinase K. They do not show any signal, since PrP^c was fully digested.

Positive samples were also all treated with proteinase K.

BSE positive bovine sample (position 2), **classical scrapie positive sample** (position 6), and **CWD positive sample** (position 4) show a typical 3 band pattern, demonstrating digestion of PrP^c and transformation of the disease specific prion protein into a proteinase resistant core fragment with reduced molecular mass following removal of the N-terminus part of the protein. The two higher bands correspond to mono and di-glycosylated forms (27-30 kDa) while the lower band corresponds to the non glycosylated form.

Ovine sample experimentally infected with BSE (position 5) is presenting a higher signal on the di-glycosylated than on the mono-glycosylated band. Nevertheless, this typical glycoprofile can not be considered as a sufficient proof of the infection of the animal with BSE. According to the Community Reference Laboratory (CRL) recommendations, a differentiation assay must be performed on this type of sample to conclude between scrapie and BSE. Please contact Bio-Rad for more information on the Bio-Rad Discriminatory test.

Atypical scrapie (e.g. Nor98) affected ovine sample (position 8) is presenting an atypical glycoprofile. A lower band is visible at approximately 12 kDa, while other upper bands are not located at the same positions compare to “typical” scrapie cases. Signal is also stronger on the lowest band than on the upper band.

Gel reading must be considered cautiously since a strong positive sample detected with the TeSeE™ WESTERN BLOT assay may hide the nearest negative or low positive sample.

Limits of the test:

A negative result means that the test sample does not contain detectable PrP^{res} by TeSeE™ WESTERN BLOT assay. However, as very low levels of PrP^{res} may not be detected, such a result does not exclude the possibility of infection.

Any sample with a negative result according to the TeSeE™ WESTERN BLOT assay interpretation criteria, following a positive rapid test result, should be tested with one of the other OIE certified confirmatory methods, Immunohistochemistry (IHC) or SAF-Immunoblot.

Any sample with a reproducible positive result according to the test interpretation criteria must be verified in accordance with current legal regulation.

8 - PRECAUTIONS

The quality of the data obtained depends on compliance with the following good laboratory practices:

- Reagents must be stored at the appropriate temperature (refer to supplier's indications).
- Do not use reagents whose shelf-life has expired.
- Do not use reconstituted proteinase K after 10 hours storage at room temperature (+18°C to +30°C).
- Do not mix or combine reagents derived from different batches of the TeSeE™ WESTERN BLOT assay during the same manipulation, with the exception of grinding tubes, reagent A, reagent B and proteinase K.
- Allow the reagents and buffers to adjust to room temperature (+18°C to +30°C) for 30 minutes before use.
- Thoroughly reconstitute reagents, avoiding any contamination.
- Do not perform the test in the presence of reactive vapors (acids, alkalines, aldehydes) or dust, which could alter the enzymatic activity of the conjugate.
- The enzymatic reaction is very sensitive to all metals or metallic ions. Consequently, no metallic element must be in contact with the conjugate.
- Only use polypropylene tubes.
- Use clean glassware, rinsed in distilled water, or preferably disposable material.
- Use a new pipette tip for each sample.
- When starting electrophoresis and transfer, check that the 2 electrodes are in contact with buffer.
- All the rinsing times must be respected to avoid any excess background noise during final staining with ECL reagent.

9 - HYGIENE AND SAFETY INSTRUCTIONS

Generally, hygiene conditions, biosafety measures and good laboratory practices must be in agreement with the recommendations of national regulatory authorities.

- All reagents of the kit are intended for use in “*in vitro*” diagnosis.
- Wear disposable gloves when handling reagents and samples and wash your hands thoroughly after handling them.
- Do not pipette by mouth.
- Use polypropylene containers to avoid broken glass.
- All the materials directly in contact with the samples and the wash solutions must be considered as contaminated.
- Avoid splashing samples or solutions containing samples.
- Contaminated surfaces must be cleaned with 20 000 ppm sodium hypochlorite solution (bleach). When the contaminating liquid is an acid, contaminated surfaces must be first neutralized with sodium hydroxide before using bleach. Surfaces must be rinsed with distilled water, dried with ethanol and wiped with absorbent paper. The material used for cleaning must be discarded in a specific container for contaminated waste.
- Samples, material and contaminated products must be eliminated after decontamination:
 - either by soaking in 1 M sodium hydroxide (final concentration) for at least 1 hour at room temperature (+18°C to +30°C),
 - or by soaking in 20 000 ppm sodium hypochlorite solution for at least 1 hour at room temperature (+18°C to +30°C),
 - or by autoclaving at 134°C minimum for at least 18 minutes, under 3 bars of pressure.

Note: never autoclave solutions containing bleach or reagent B.

- All operations involved in Transmissible Spongiform Encephalopathy (TSE) screening tests are subject to regulations and must be performed in an isolated, limited and controlled access laboratory devoted exclusively to this activity. A laboratory coat, overshoes, gloves, mask with visor or simple mask with safety glasses are required to ensure the operator’s safety.
- Operators must receive specific training concerning the risks related to TSEs agents or prions and the validated modes of decontamination for unconventional agents. Biosafety measures must be in agreement with recommendations of regular authorities of the country.

- Neutralize and/or autoclave all wash solutions or wash wastes or any liquid containing biological samples prior to their elimination.
- Reagent B is a dangerous substance classified as nocive (> 25% alcohol) according to European legislation.
- Reagents containing 0.1% ProClin® 300 are classified as irritating preparations according to European legislation.

Xn
(Alcohol > 25%)
(0.1% ProClin® 300)

R : 10-22-37/38-41-43-67 Flammable. Harmful if swallowed. Irritating to respiratory system and skin. Risk of serious damage to eyes. May cause sensitisation by skin contact. Inhalation of vapour may cause drowsiness and dizziness.

S : 7/9-13-26-28-37/39-46 Keep container tightly closed and in a well ventilated place. Keep away from food, drink and animal feed. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin wash immediately with plenty of water. Wear suitable protecting clothings, gloves and eye/face protection. If swallowed, seek medical advice immediately and show this container or label.

10 - REFERENCES

1. S.B. PRUSINER (1991)
Molecular biology of prion diseases - *Science* 252:1515-1522.
2. J.B. KATZ, J.G. PEDERSEN, A.L. JENNY and W.D. TAYLOR (1992)
Assessment of Western immunoblotting for the confirmatory diagnosis of ovine scrapie and bovine spongiform encephalopathy (BSE) - *Journal of Veterinary Diagnostic Investigations*, 4, 447-449.
3. G.A.H. WELLS, Y.I. SPENCER and M. HARITANI (1994)
Configuration and topographic distribution of PrP in the central nervous system in bovine spongiform encephalopathy: an immunohistochemical study. In: *Slow Infections of the Central Nervous System*, J. BJORNSSON, R.I. CARP, A. LOVE and M. WISNIEWSKI - Eds *The New York Academy of Sciences*, pp 350-352.
4. J.P. DESLYS, E. COMOY, S. HAWKINS, S. SIMON, H. SCHIMMEL, G. WELLS, J. GRASSI and J. MOYNAGH (2001)
Screening slaughtered cattle for BSE - *Nature*: 409; 476-477.
5. S. BENESTAD, P. SARRADIN, B. THU, J. SCHÖNHEIT, M. TRANULIS and B. BRATBERG (2003)
Cases of scrapie with unusual features in Norway and designation of new type, Nor98. *Veterinary Record* 153, 202-208.
6. A. BUSCHMANN, G. LÜHKEN, J. SCHULTZ, G. ERHARDT and M.-H. GROSCHUP (2004)
Neuronal accumulation of abnormal prion protein in sheep carrying a scrapie-resistant genotype (PrP^{ARRV/ARRV}). *Journal of General Virology* 85, 2727-2733.
7. L. ORGE, A.GALO, C.MACHADO, C. LIMA, C. OCHOA, J. SILVA, M. RAMOS and J.-P. SIMAS (2004)
Identification of putative atypical scrapie in sheep in Portugal. *Journal of General Virology* 85, 3487-3491.
8. A. BUSCHMANN, A.-G. BIACABE, U. ZIEGLER, A. BENCSIK, J.-Y. MADEC, G. ERHARDT, G. LÜHKEN, T. BARON and M.-H. GROSCHUP (2004)
Atypical scrapie cases in Germany and France are identified by discrepant reaction patterns in BSE rapid tests. *Journal of Virological Methods* 117, 27-36.

9. H. DE BOSSCHERE, S. ROELS, S. BENESTAD, E. VANOPDENBOSCH (2004). Scrapie case similar to Nor98 diagnosed in Belgium via active surveillance. *Veterinary Record* 155, 707-708.
10. H. ONNASCH, H. M. GUNN, B.J. BRADSHAW, S. BENESTAD, H.F. BASSETT (2004). Two Irish cases of scrapie resembling Nor98. *Veterinary Record* 155, 636-637.
11. S. BENESTAD, P. SARRADIN, J.-M. BILHEUDE, J. GRASSI, H. LAUDE, O. ANDREOLETTI, T. MOLDAL and B. BRATBERG. Are there gold standard methods for the diagnosis of TSE ? Nor98 scrapie cases: atypical cases and their challenging diagnosis. Second International Chronic Wasting Disease Symposium, Madison - Wisconsin, USA.
12. BIACABE, A-G., LAPLANCHE, J-L., RYDER, S. AND BARON, T. (2004). Distinct molecular phenotypes in bovine prion diseases, *EMBO Reports* 5, 110-114.
13. CASALONE, C., ZANUSSO, G., ACUTIS, P., FERRARI, S., CAPUCCI, I., TAGLIAVINI, F., MONACO, S. and CARAMELLI, M. (2004). Identification of a second bovine amyloidotic spongiform encephalopathy: Molecular similarities with sporadic Creutzfeldt-Jakob disease. *PNAS* 101, 3065-70.
14. BUSCHMANN A; GRETZSCHEL A; BIACABE A-G; SCHIEBEL K; CORONA C; HOFFMANN C; EIDEN M; BARON T; CASALONE C; GROSCHEP M H (2006). Atypical BSE in Germany-proof of transmissibility and biochemical characterization. *Veterinary microbiology* 2006;117(2-4):103-16.
15. ARSAC, J.-N., ANDREOLETTI, O., BILHEUDE, J.-M., LACROUX, C., BENESTAD, S. L., AND BAARON, T. (2007). Similar biochemical signatures and prion protein genotypes in atypical scrapie and Nor98 cases, France and Norway. *Emerging Infectious Diseases* 13(1), 58-65.
16. E. URO-COSTE, H CASSARD, S. SIMON, S. LUGAN, J.-M. BILHEUDE, A. PERRET-LIAUDET, J.-W. IRONSIDE, S. HAIK, C. BASSET-LEOBON, C. LACROUX, K. PEOCH, N. STREICHENBERGER, J. LANGEVELD, M.-W. HEAD, J. GRASSI, J.-J. HAUW, F. SCHELCHER, M.-B. DELISLE, O. ANDREOLETTI (2008) Beyond PrP^{res} Type 1/Type 2 Dichotomy in Creutzfeldt-Jakob disease. *Plos Pathogens* volume 4, issue 2, e1000029.

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10 - BIBLIOGRAPHIE

1. S.B. PRUSINER (1991)
Molecular biology of prion diseases - *Science* 252:1515-1522.
2. J.B. KATZ, J.G. PEDERSEN, A.L. JENNY and W.D. TAYLOR (1992) Assessment of Western immunoblotting for the confirmatory diagnosis of ovine scrapie and bovine spongiform encephalopathy (BSE) - *Journal of Veterinary Diagnostic Investigations*, 4, 447-449.
3. G.A.H. WELLS, Y.I. SPENCER and M. HARITANI (1994)
Configuration and topographic distribution of PrP in the central nervous system in bovine spongiform encephalopathy: an immunohistochemical study. In: *Slow Infections of the Central Nervous System*, J. BJORNSSON, R.I. CARP, A. LOVE and M. WISNIEWSKI - Eds *The New York Academy of Sciences*, pp 350-352.
4. J.P. DESLYS, E. COMOY, S. HAWKINS, S. SIMON, H. SCHIMMEL, G. WELLS, J. GRASSI and J. MOYNAGH (2001)
Screening slaughtered cattle for BSE - *Nature*: 409; 476-477.
5. S. BENESTAD, P. SARRADIN, B. THU, J. SCHÖNHEIT, M. TRANULIS and B. BRATBERG (2003)
Cases of scrapie with unusual features in Norway and designation of new type, Nor98. *Veterinary Record* 153, 202-208.
6. A. BUSCHMANN, G. LÜHKEN, J. SCHULTZ, G. ERHARDT and M.-H. GROSCHUP (2004)
Neuronal accumulation of abnormal prion protein in sheep carrying a scrapie-resistant genotype (PrP^{ARR/ARR}). *Journal of General Virology* 85, 2727-2733.
7. L. ORGE, A.GALO, C.MACHADO, C. LIMA, C. OCHOA, J. SILVA, M. RAMOS and J.-P. SIMAS (2004)
Identification of putative atypical scrapie in sheep in Portugal. *Journal of General Virology* 85, 3487-3491.
8. A. BUSCHMANN, A.-G. BIACABE, U. ZIEGLER, A. BENCSIK, J.-Y. MADEC, G. ERHARDT, G. LÜHKEN, T. BARON and M.-H. GROSCHUP (2004)
Atypical scrapie cases in Germany and France are identified by discrepant reaction patterns in BSE rapid tests. *Journal of Virological Methods* 117, 27-36.

9. H. DE BOSSCHERE, S. ROELS, S. BENESTAD, E. VANOPDENBOSCH (2004). Scrapie case similar to Nor98 diagnosed in Belgium via active surveillance. *Veterinary Record* 155, 707-708.
10. H. ONNASCH, H. M. GUNN, B.J. BRADSHAW, S. BENESTAD, H.F. BASSETT (2004). Two Irish cases of scrapie resembling Nor98. *Veterinary Record* 155, 636-637.
11. S. BENESTAD, P. SARRADIN, J.-M. BILHEUDE, J. GRASSI, H. LAUDE, O. ANDREOLETTI, T. MOLDAL and B. BRATBERG. Are there gold standard methods for the diagnosis of TSE ? Nor98 scrapie cases: atypical cases and their challenging diagnosis. Second International Chronic Wasting Disease Symposium, Madison - Wisconsin, USA.
12. BIACABE, A-G., LAPLANCHE, J-L., RYDER, S. AND BARON, T. (2004). Distinct molecular phenotypes in bovine prion diseases, *EMBO Reports* 5, 110-114.
13. CASALONE, C., ZANUSSO, G., ACUTIS, P., FERRARI, S., CAPUCCI, I., TAGLIAVINI, F., MONACO, S. and CARAMELLI, M. (2004). Identification of a second bovine amyloidotic spongiform encephalopathy: Molecular similarities with sporadic Creutzfeldt-Jakob disease. *PNAS* 101, 3065-70.
14. BUSCHMANN A; GRETZSCHEL A; BIACABE A-G; SCHIEBEL K; CORONA C; HOFFMANN C; EIDEN M; BARON T; CASALONE C; GROSCHUP M H (2006). Atypical BSE in Germany-proof of transmissibility and biochemical characterization. *Veterinary microbiology* 2006;117(2-4):103-16.
15. ARSAC, J.-N., ANDREOLETTI, O., BILHEUDE, J.-M., LACROUX, C., BENESTAD, S. L., AND BAARON, T. (2007). Similar biochemical signatures and prion protein genotypes in atypical scrapie and Nor98 cases, France and Norway. *Emerging Infectious Diseases* 13(1), 58-65.
16. E. URO-COSTE, H CASSARD, S. SIMON, S. LUGAN, J.-M. BILHEUDE, A. PERRET-LIAUDET, J.-W. IRONSIDE, S. HAIK, C. BASSET-LEOBON, C. LACROUX, K. PEOCH, N. STREICHENBERGER, J. LANGEVELD, M.-W. HEAD, J. GRASSI, J.-J.HAUW, F. SCHELCHER, M.-B. DELISLE, O. ANDREOLETTI (2008) Beyond PrP^{res} Type 1/ Type 2 Dichotomy in Creutzfeldt-Jakob disease. *Plos Pathogens* volume 4, issue 2, e1000029.

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1. S.B. PRUSINER (1991)
Molecular biology of prion diseases - *Science* 252:1515-1522.
2. J.B. KATZ, J.G. PEDERSEN, A.L. JENNY and W.D. TAYLOR (1992) Assessment of Western immunoblotting for the confirmatory diagnosis of ovine scrapie and bovine spongiform encephalopathy (BSE) - *Journal of Veterinary Diagnostic Investigations*, 4, 447-449.
3. G.A.H. WELLS, Y.I. SPENCER and M. HARITANI (1994) Configuration and topographic distribution of PrP in the central nervous system in bovine spongiform encephalopathy: an immunohistochemical study. In: *Slow Infections of the Central Nervous System*, J. BJORNSSON, R.I. CARP, A. LOVE and M. WISNIEWSKI - Eds *The New York Academy of Sciences*, pp 350-352.
4. J.P. DESLYS, E. COMOY, S. HAWKINS, S. SIMON, H. SCHIMMEL, G. WELLS, J. GRASSI and J. MOYNAGH (2001) Screening slaughtered cattle for BSE - *Nature*: 409; 476-477.
5. S. BENESTAD, P. SARRADIN, B. THU, J. SCHÖNHEIT, M. TRANULIS and B. BRATBERG (2003) Cases of scrapie with unusual features in Norway and designation of new type, Nor98. *Veterinary Record* 153, 202-208.
6. A. BUSCHMANN, G. LÜHKEN, J. SCHULTZ, G. ERHARDT and M.-H. GROSCHUP (2004) Neuronal accumulation of abnormal prion protein in sheep carrying a scrapie-resistant genotype (PrP^{PARR/ARR}). *Journal of General Virology* 85, 2727-2733.
7. L. ORGE, A.GALO, C.MACHADO, C. LIMA, C. OCHOA, J. SILVA, M. RAMOS and J.-P. SIMAS (2004) Identification of putative atypical scrapie in sheep in Portugal. *Journal of General Virology* 85, 3487-3491.
8. A. BUSCHMANN, A.-G. BIACABE, U. ZIEGLER, A. BENCSIK, J.-Y. MADEC, G. ERHARDT, G. LÜHKEN, T. BARON and M.-H. GROSCHUP (2004) Atypical scrapie cases in Germany and France are identified by discrepant reaction patterns in BSE rapid tests. *Journal of Virological Methods* 117, 27-36.

9. H. DE BOSSCHERE, S. ROELS, S. BENESTAD, E. VANOPDENBOSCH (2004). Scrapie case similar to Nor98 diagnosed in Belgium via active surveillance. *Veterinary Record* 155, 707-708.
10. H. ONNASCH, H. M. GUNN, B.J. BRADSHAW, S. BENESTAD, H.F. BASSETT (2004). Two Irish cases of scrapie resembling Nor98. *Veterinary Record* 155, 636-637.
11. S. BENESTAD, P. SARRADIN, J.-M. BILHEUDE, J. GRASSI, H. LAUDE, O. ANDREOLETTI, T. MOLDAL and B. BRATBERG. Are there gold standard methods for the diagnosis of TSE ? Nor98 scrapie cases: atypical cases and their challenging diagnosis. Second International Chronic Wasting Disease Symposium, Madison - Wisconsin, USA.
12. BIACABE, A-G., LAPLANCHE, J-L., RYDER, S. AND BARON, T. (2004). Distinct molecular phenotypes in bovine prion diseases, *EMBO Reports* 5, 110-114.
13. CASALONE, C., ZANUSSO, G., ACUTIS, P., FERRARI, S., CAPUCCI, I., TAGLIAVINI, F., MONACO, S. and CARAMELLI, M. (2004). Identification of a second bovine amyloidotic spongiform encephalopathy: Molecular similarities with sporadic Creutzfeldt-Jakob disease. *PNAS* 101, 3065-70.
14. BUSCHMANN A; GRETZSCHEL A; BIACABE A-G; SCHIEBEL K; CORONA C; HOFFMANN C; EIDEN M; BARON T; CASALONE C; GROSCHUP M H (2006). Atypical BSE in Germany-proof of transmissibility and biochemical characterization. *Veterinary microbiology* 2006;117(2-4):103-16.
15. ARSAC, J.-N., ANDREOLETTI, O., BILHEUDE, J.-M., LACROUX, C., BENESTAD, S. L., AND BAARON, T. (2007). Similar biochemical signatures and prion protein genotypes in atypical scrapie and Nor98 cases, France and Norway. *Emerging Infectious Diseases* 13(1), 58-65.
16. E. URO-COSTE, H CASSARD, S. SIMON, S. LUGAN, J.-M. BILHEUDE, A. PERRET-LIAUDET, J.-W. IRONSIDE, S. HAIK, C. BASSET-LEOBON, C. LACROUX, K. PEOCH, N. STREICHENBERGER, J. LANGEVELD, M.-W. HEAD, J. GRASSI, J.-J.HAUW, F. SCHELCHER, M.-B. DELISLE, O. ANDREOLETTI (2008) Beyond PrP^{res} Type 1/ Type 2 Dichotomy in Creutzfeldt-Jakob disease. *Plos Pathogens* volume 4, issue 2, e1000029.

APPENDIX

OIE Procedure for Registration of Diagnostic Kits Abstract sheet

Name of the diagnostic kit: TeSeE™ WESTERN BLOT

Manufacturer: Bio-Rad

OIE Approval number: 20090105

Date of Registration: May 2009

Disease: Prion associated diseases (Scrapie, Bovine Spongiform Encephalopathy, Chronical Waste Disease)

Pathogen Agent: Abnormal prion protein PrPRes

Type of Assay: The TeSeE™ WESTERN BLOT kit is based on the Western blotting technique and allows detection of abnormal prion protein PrPRes after protease digestion.

Purpose of Assay: Certified by the OIE in May 2009 as fit for the post-mortem detection of transmissible spongiform encephalopathies (TSEs) in cattle (bovine spongiform encephalopathy, BSE), in ovines and caprines (BSE and scrapie), and in cervids (Chronic Wasting Disease, CWD), and for the following purposes:

1. To confirm TSE suspected positive samples detected at the screening laboratories in countries with active/passive surveillance programmes. Any sample with a negative result according to the TeSeE™ WESTERN BLOT assay interpretation criteria, following a positive rapid test result, should be tested with one of the other OIE certified confirmatory methods, Immunohistochemistry (IHC) or Scrapie-associated fibrils (SAF)-Immunoblot;
2. To confirm the prevalence of infection with one of the TSE associated diseases (BSE, scrapie, CWD) in the context of an epidemiological survey in a low prevalence country;
3. To estimate prevalence of infection to facilitate risk analysis (e.g. surveys, implementation of disease control measures) and to assist the demonstration of the efficiency of eradication policies.

Species and Specimen: Bovine (obex), Cervids (obex and lymph nodes), and ovine and caprine (obex)

1. Information on the kit

Information is available on Bio-Rad Website, at www.bio-rad.com

2. Summary of validation studies

ANALYTICAL CHARACTERISTICS

Calibration:

Since no reference standard is available, the TeSeE™ WESTERN BLOT assay was calibrated against the Bio-Rad TeSeE rapid test (TeSeE immunoassay: sandwich ELISA format), on a panel of Immunohistochemistry (IHC) confirmed positive samples, including weak positive and atypical scrapie samples. The TeSeE™ WESTERN BLOT could detect all positive samples.

Repeatability:

Both within run and between run repeatability were evaluated.

Within run: The results were perfectly reproducible when testing 32 times within the same day, a panel of 1 negative and 3 positive samples. Within run repeatability was confirmed by the Reference Laboratory for European Union (CRL), VLA, Weybridge, UK, within assay study.

Between run: The results were perfectly reproducible when testing a panel of 1 negative and 3 positive samples, in duplicate, for 20 days.

Between serial repeatability is verified for every reagent batch produced through the final quality control procedure performed in the manufacturing site with a specifically designed Quality Control sample panel.

Analytical specificity:

As specimens collected on animals infected both prion and other pathogen are difficult to obtain, analytical specificity was not studied.

Analytical sensitivity:

The TeSeE™ WESTERN BLOT assay was demonstrating the same sensitivity than the Bio-Rad TeSeE S/G Rapid test (ELISA) and than the OIE scrapie-associated fibrils (SAF) blot method.

DIAGNOSTIC CHARACTERISTICS:

Field samples from passive or active surveillance programs were used for the evaluation studies of the TeSeE™ WESTERN BLOT assay. The status (positive or negative) of these samples had been determined with one of the European Commission approved screening methods and/or IHC reference method. The ability of the TeSeE™ WESTERN BLOT assay to confirm their status was evaluated. The same evaluation approach was used in all evaluations performed.

- Study I Internal evaluation at Bio-Rad, France
- Analytical sensitivity on ovine and bovine characterized samples from veterinary reference laboratories
 - Diagnostic specificity on ovine and bovine samples from French slaughterhouses (animal without any clinical sign) classified as negative with a rapid screening tests (TeSeE™ ELISA)
 - Diagnostic sensitivity on ovine samples from French flocks with natural cases of TSE tested with the TeSeE™ ELISA and IHC;
- Study II External evaluation at CRL, VLA, Weybridge, UK
- Diagnostic specificity and sensitivity on bovine samples from UK passive surveillance
 - Diagnostic specificity and sensitivity retrospective study on bovine samples from UK
 - Diagnostic specificity and sensitivity on ovine samples from UK active and passive surveillance;
- Study III External study at CODA-CERVA, Belgium: prospective study on bovine samples from routine screening and retrospective study on previously BSE positive samples from the laboratory routine;
- Study IV External study at IZW Berlin, Germany: retrospective study on 104 cervid samples (deer and elk) from Germany, USA and Canada, characterized by IHC including 64 CWD positive samples;
- Study V External study at AFSSA, France: retrospective study on scrapie positive ovine samples chosen randomly among the scrapie diagnosed cases between 2002 and 2004 in France (positive with the French reference AFSSA Western blot method and with the TeSeE™ Sheep/Goat ELISA);
- Study VI External retrospective study at AFSSA, France: Arzac et al, Acta Neuropathol., July 2007: retrospective study on BSE positive obex collected in France during active and passive surveillance between 2001 and 2003 (confirmed positive with the AFSSA Western blot and analysed with IHC in case of discrepancy between rapid test and Western blot);
- Study VII External retrospective study at AFSSA, France: Arzac et al, EID January 2007: retrospective study on 54 atypical scrapie positive samples detected by the rapid test TeSeE™ ELISA, collected in France during the scrapie surveillance program between 2002 and 2004;

- Study VIII External retrospective study on BSE positive samples (including atypical BSE cases) at CRL, VLA, Weybridge: Terry et Al., Veterinary Record, June 2007: retrospective study on 5 BSE positive confirmed samples (including one French H type from AFSSA and one H type from UK) and 1 scrapie positive confirmed sample;
- Study IX CFIA Ottawa (external study): retrospective study on 50 positive elk brain samples, 40 positive White Tailed Deer brain, 65 positive White Tailed Deer retropharyngeal lymph node samples;
- Study X USDA Ames Iowa (external study): retrospective study on 53 positive deer retropharyngeal lymph node samples.

THRESHOLD DETERMINATION

The threshold was determined according to the following criteria:

Negative samples: They do not show any signal, since PrPc was totally digested by proteinase K.

Positive samples:

- Classical positive samples: They show a typical 3 band pattern, demonstrating digestion of PrPc and transformation of the disease specific prion protein into a proteinase resistant core fragment with reduced molecular mass following removal of N-terminus part of the protein. The two higher bands correspond to mono and di-glycosylated forms (27-30 kDa) while the lower band corresponds to the un-glycosylated form of the protein.
- Nor 98 positive samples: They show an atypical glycoprofile. A lower band is visible at approximately 12kDa, while other upper bands are not located at the same positions compare to classical scrapie cases.
- Other samples that should be repeated: The samples for which only the di-glycosylated band is detected and the samples for which only the di-glycosylated band and the mono-glycosylated band are detected.

Diagnostic sensitivity (DSn) and specificity (DSp) (with 95% confidence limits [CI])

TeSeE™ WESTERN BLOT versus VLA Hybrid WB		Bovine
Diagnostic sensitivity	N	223
	DSn	100 %
	CI	98,4 %–100 %
Diagnostic specificity	N	50
	DSp	100 %
	CI	92,9 %–100 %

TeSeE™ WESTERN BLOT versus AFSSA WB		Bovine
Diagnostic sensitivity	N	52
	DSn	100 %
	CI	93,1 %–100 %
Diagnostic specificity	N	25
	DSp	84 %
	CI	63,9 %–95,5 %

TeSeE™ WESTERN BLOT versus IHC		Ovine
Diagnostic sensitivity	N	211
	DSn	99 %
	CI	96,4 % - 99,9 %
Diagnostic specificity	N	71
	DSp	100 %
	CI	94,9 %–100 %

TeSeE™ WESTERN BLOT versus SAF Immunoblot OIE		Ovine
Diagnostic sensitivity	N	73
	DSn	100 %
	CI	95,1 %–100 %
Diagnostic specificity	N	61
	DSp	100 %
	CI	94,1 %–100 %

TeSeE™ WESTERN BLOT versus SAF Immunoblot OIE		Ovine
Diagnostic sensitivity	N	73
	DSn	100 %
	CI	95,1 %–100 %
Diagnostic specificity	N	61
	DSp	100 %
	CI	94,1 %–100 %

(*): the 54 discrepant samples were all confirmed positive by IHC (atypical scrapie)

TeSeE™ WESTERN BLOT versus AFSSA WB		Ovine
Diagnostic sensitivity	N DSn CI	40 100 % 91,2 %–100 %
Diagnostic specificity	N DSp CI	0 NC NC

TeSeE™ WESTERN BLOT versus IHC		Cervids
Diagnostic sensitivity	N DSn CI	272 100 % 98.65 %–100 %
Diagnostic specificity	N DSp CI	40 100 % 91,2 %–100 %

REPRODUCTIBILITY

The CRL regularly issues proficiency testing panels for a range of TSE tests. The samples consist of homogenised brain material presented as 50% tissue/water homogenates. The summarised results for confirmatory blotting rounds issued in 2007 are given below.

- CRL Scrapie confirmatory blotting Proficiency test round

A panel of 5 samples (1 strong positive, 2 weak positives [all classical scrapie] and 2 negatives) was issued to 23 European laboratories in November 2007.

10 sets of results were returned by laboratories using the BioRad TeSeE™ WESTERN BLOT. All of these 10 laboratories identified all samples correctly.

- CRL BSE confirmatory blotting Proficiency test round

A panel of 5 samples (1 strong positive, 1 medium positive, 1 weak positive [all classical BSE] and 2 negatives) was issued to 21 European laboratories in October 2007.

Eight sets of results were returned by laboratories using the BioRad TeSeE™ WESTERN BLOT. All of these 8 labs identified all samples correctly.

In addition, reproducibility of the TeSeE™ WESTERN BLOT was evaluated with the Bio-Rad control panel including 2 negative ovine samples, 1 negative bovine sample, 1 CWD positive sample, 3 scrapie positive samples and 1 BSE positive sample. The panel was provided to 3 TSE Reference laboratories: VLA in UK, CFIA in Canada and NVI in Norway. The results from the 3 sites were compared to the results obtained in the Bio-Rad facilities. All of these 3

laboratories identified all samples correctly (the BSE positive sample was not evaluated by the CFIA ACIA laboratory in Ottawa since this laboratory is not a Reference laboratory for BSE in Canada).

APPLICATIONS

The TeSeE™ WESTERN BLOT assay is only for use in National Reference Laboratories (NRLs) and Community Reference laboratories (CRLs). The TeSeE™ WESTERN BLOT assay can be used for different types of applications:

- **as a confirmatory assay** for the confirmation of TSE suspected cases, in cattle (BSE), small ruminants e.g. sheep/goat (BSE/Scrapie) and cervids (CWD). The diagnostic regimen includes several steps at both screening laboratories and Reference laboratories.
- at the screening laboratory: initial screening with one of the approved rapid test. If the sample turns positive (= initial reactive), the test needs to be repeated in duplicate with the same assay and with the same sample (sample homogenate). If one of the two duplicates or the two turn positive, the sample is then considered as “suspect”. All remaining material (sample homogenate + remaining tissues) are sent to the Reference laboratory for confirmation.
- at the reference laboratory: A negative result with the TeSeE™ WESTERN BLOT means that the test sample does not contain detectable PrPres by TeSeE™ WESTERN BLOT assay. However, as very low levels of PrPres may not be detected, such a result does not exclude the possibility of infection.

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- as a rapid method for the determination of the prevalence of infection with one of the TSE associated diseases (BSE, scrapie, CWD) in the context of an epidemiological survey in a low prevalence country.

REFERENCES

1. J.N. Arsac, A.G Biacabe, J.Nicollo, A.Bencsik and T. Baron. Biochemical identification of bovine spongiform encephalopathies in cattle. Acta Neuropathologica 2007, vol 114:509-516
2. J.N. Arsac, O. Andreoletti, J.M. Bilheude, C. Lacroux, S. L. Benestad and T. Baron. Similar Biochemical Signatures and Prion Protein Genotypes in Atypical Scrapie and Nor98 Cases, France and Norway. Emerging Infectious Diseases –Vol. 13, No. 1, January 2007, 58-65.
3. Terry, L.A., Jenkins, R., Thorne, L., Everest, S.J., Chaplin, M.C., Davis, L.A. And Stack, M.J. (2007) First Case Of H-Type Bovine Spongiform Encephalopathy Identified In Great Britain. Veterinary Record 160: 873 – 874.
4. Biacabe A.G, Morignat E, Vulin J, Calavas D, Baron TG.M; February 2008 Atypical Bovine Spongiform Encephalopathies, France, 2001-2007. Emerging Infectious Diseases Vol.14, n°2

Purchase of Criterion XT Bis-Tris gels and XT MOPS running buffer is accompanied by a limited license under U.S. Patent Numbers 6,143,154; 6,096,182; 6,059,948; 5,578,180; 5,922,185; 6,162,338; and 6,783,651, and corresponding foreign patents.

ECL is a trademark of GE Healthcare. Eppendorf is a trademark of Eppendorf AG. Kodak is a trademark of Eastman Kodak Company.

Bio-Rad

3, boulevard Raymond Poincaré
92430 MARNES-LA-COQUETTE - FRANCE
Tel.: +33 1 47 95 60 00
Fax.: +33 1 47 41 91 33

14-1032 0614
Code : 862231